

# Programmable and Stimuli-Responsive Hydrogels: AI-Assisted Design for Advanced Drug Delivery and Tissue Engineering

Foram Patel\*, Hitesh Jain, D. B. Meshram

*Pioneer Pharmacy College, Sayajipura, Vadodara*

*\*Corresponding Author*

Date of Submission: 11-05-2026

Date of Acceptance: 25-05-2026

## Abstract

The recent introduction of hydrogels capable of being programmed for active response to environmental stimuli has been a game-changer in terms of biomedical applications. AI-guided design, joined with multiscale simulations, helps researchers to get rid of the trial-and-error approach of traditional development methods, resulting in hydrogels that are good for targeted, site-specific drug delivery and advanced 4D bioprinting tissue engineering.

**Keyword:** Programmable Hydrogels, Stimuli-Responsive Hydrogels, Artificial Intelligence

## I. Introduction

Smart/Stimuli Responsive Hydrogels (Programmable Hydrogels) have a reactive (Dynamic) Closed-Loop Response, versus the Fixed (Static), linear Single Response of Conventional Hydrogels. CV Hydrogel response characteristics: Hydrogels have fixed water absorption capacity (swelling) regardless of outside forces acting upon them; Hydrogels cannot control the degree they transition between gel and solution; Hydrogels possess lower permeability than Smart Hydrogel systems; Hydrogels (Standalone) Stronger than Smart Hydrogel systems; Smart Hydrogels respond only to external stimuli on demand by design to provide adaptive (multiple) form and function (Smart Hydrogels can provide multiple (adaptive) forms and functions). An important difference between Programmable Hydrogels vs. Conventional Hydrogel systems is the capacity for Programmable Hydrogel systems to run inside a Closed-Loop Feedback Control System: Programmable Hydrogel systems are usually able to make sense of something in the environment (environmental or biological signal) and begin a specific Bio-Function, for example, release an agent to treat a disease, based on the response to the signals received.<sup>(1)</sup>

Increasingly significant in the academic realm, AI-assisted development of hydrogels will use ML-based methodologies for creating and improving hydrogel properties. AI is enhancing the production of responsive intelligent hydrogel. The application of computational modelling activates AI to make available a scientific basis for the development of hydrogels based on structure-property correlations; as a result, it allows the prediction of hydrogels' behaviour under many conditions. AI hydrogel materials are intelligent and can assist in developing multimodal cancer therapy by providing improved monitoring capabilities for the tumour microenvironment (TME) as well as a means of dynamically modifying it. With decreased time and cost associated with R&D facilitated through the use of AI technology applications, the use of AI will enable hydrogels to be produced in an automatable manner for use in transitioning from laboratory settings to clinical settings.<sup>(2,3,4,5,7)</sup>

## Types of Programmable Hydrogels

- **Temperature Programmable Hydrogels**

Temperature Programmable hydrogels change their properties and configurations with temperature changes, which makes a difference in how their polymers are arranged and connected.<sup>(6)</sup>

- **Mechanically Programmable Hydrogels**

Mechanically Programmable hydrogels allow for physical input to change their behaviour based on input from mechanical sources (like light or force) by converting these external stimuli into physicochemical energy so that they can reform and reprogram their structure.<sup>(6)</sup>

- **Shape Programmable Hydrogels**

Shape programmable hydrogels are smart soft materials that show so many complex deformations by responding to environmental stimuli, for example, a change like temperature, pH, introduction of light, magnetism, electrical current, or the addition of water.<sup>(6)</sup>

- **Gene Programmable Hydrogels**

Through the complementary hybridisation of nucleic acids, gene-programmable hydrogels are able to react to particular genetic instructions or sequences.<sup>(6)</sup>

- **Colour Programmable Hydrogels**

Colour Programmable hydrogels are biomimetic, and it is used as camouflage or sensing materials. Still, conventional colour-changing hydrogel typically has poor correctness due to a lack of confidence around the occurrence of colour changes, as these chemical processes are responsible.<sup>(6)</sup>

- **Degradable Programmable Hydrogels**

Hydrogels utilise numerous stimulants such as light, pH, enzymes, and CRISPR systems to yield programmable degradation rates and structural breakdown profiles, which allow them to be safely and precisely used in biomedical applications.<sup>(6)</sup>

- **Adhesive Programmable Hydrogels**

Traditional hydrogels have so many limitations that is why they compromise their capability to cause biological integrity and provide suitable adhesion to wearables, containing poor wet adhesion and poor control of the material substrate.<sup>(6)</sup>

- **Function/Response Programmable Hydrogels**

Function Programmable hydrogels are an outstanding choice for biomaterial development because they achieve an ideal balance of structural durability, mechanical performance and dynamic reactivity by showing highly adjustable properties.<sup>(6)</sup>

- **Transparent Programmable Hydrogels**

Programmable hydrogels are able to respond to multiple environmental stimuli such as light, temperature, and pH in a manner that changes their structural characteristics, transparency, or colour, so that they can be used to create advanced, multilayer anti-counterfeiting encryption that is far more sophisticated than traditional single-stimulus methods.<sup>(6)</sup>

### Types of stimuli – Responsive Hydrogels

- **Temperature-Responsive Hydrogels**

Temperature-Responsive hydrogels are those that exhibit a marked change in physical state at a threshold temperature called a critical solution temperature, thus categorised as lower critical solution temperatures (LCST) or upper critical solution temperatures (UCST).<sup>(7)</sup>

- **pH – Responsive Hydrogels**

The pH-Responsive hydrogel conforms by swelling, contracting or dissolving when the pH of the environment around the hydrogel changes due to the presence of ionizable functional groups, such as

hydrazones, imines, esters and amides in the polymer chains.<sup>(7)</sup>

- **Enzyme – Responsive Hydrogels**

As cancer progresses, the production of matrix metalloproteinases (MMPs) is upregulated by a considerable amount; MMPs are zinc-dependent proteases that specifically hydrolyse the extracellular matrix component of the tumour and cause the ECM to be degraded, and disrupt the adhesion of cells, which in turn promotes tumour invasion and metastasis.<sup>(7)</sup>

- **Light – Responsive Hydrogels**

Photo-Responsive hydrogels have photo-sensitive materials with differing sensitivities to different wavelengths of light; for that reason, they need high-energy UV light to go through roughly one micrometre.<sup>(7)</sup>

- **Electric Field-Responsive Hydrogels**

Electric-field-responsive hydrogels contain charged groups or add conductive or magneto-sensitive nanoparticles to make the hydrogel polymer matrix capable of experiencing ionic migration, electrochemical reactions, swelling, shrinking, & deforming when an electric field is applied.<sup>(1)</sup>

- **Magnetic Field-Responsive Hydrogels**

The showing of magnetic nanoparticles in a magnetic field-responsive hydrogel causes mechanical movement and heating of the hydrogel as a result of submission to a magnetic field.<sup>(1)</sup>

- **Glucose-Responsive Hydrogels**

The presence of glucose at physiological pH leads to the alteration of the permeability & degradation of glucose-sensitive hydrogels via enzymatic degradation, by using phenylboronic acid (PBA), or via concanavalin A (ConA)-mediated degradation.<sup>(1)</sup>

### Characterisation methods of smart responsive hydrogels

- **Swelling**

Mass swelling ratio, water content, and volume swelling ratio are used for physical evaluation of the swelling properties of hydrogels.<sup>(7)</sup>

- **Macroscopic Mechanical Properties**

Hydrogels are subjected to compression tests for their ability to support low-stress loads (that is, their ability not to expel water) while tensile tests are used to determine their toughness, strength, and elongation (where double-network hydrogels can stretch as much as 1000%).<sup>(7)</sup>

- **Rheology**

Dynamic Rheology is a better method of characterising the viscoelasticity of hydrogels and for studying their in-situ sol-to-gel transition kinetics.<sup>(7)</sup>

• **Chemical structure characterisation**

Hydrogel structures are typically evaluated by using a combination of spectroscopic and scattering techniques. FTIR recognize the existence of functional groups and evaluates the degree of cross-linking; NMR give extra information regarding the uniformity of the hydrogel network, conversion rates of the network from sol to gel and the density of cross-linking between the polymer chains.<sup>(7)</sup>

**AI-Driven Innovation in Programmable Hydrogel**

Material composition, formulation and environmental stimuli interplay, which is why the growth of the hydrogel was slow. To find a way around this issue, machine learning (ML) and artificial intelligence (AI) simulation have transitioned

research from being based on a trial-and-error approach to a data-driven precision methodology by pattern recognition and the optimisation of formulations with high loyalty and quality in comparison with dynamic, macroscopic responses. As an extra, AI-driven multiscale simulations that bridge molecular dynamics and continuum mechanics have been used to gain detailed awareness into the structure and function relationships of the hydrogels. These two practicalities have enabled the rational design and development of high-performance hydrogels, consequently making their transition into modern biomedical applications, for example, bioprinting, targeted drug delivery, and precision oncology.<sup>(7,8,9,10)</sup>

**Benefits and Limitations of Smart Hydrogels in Drug Delivery**

• **Benefits**

Key Features	Detailed Description & Therapeutics Impact
Controlled and Sustained Release	Smart hydrogel gives controlled, targeted delivery of medicine based on external variables, for example, temperature and pH, while reducing the frequency with which medication must be administered.
Localised Drug Delivery	Hydrogel can be used for topical application for local management of inflamed tissue.
Protection of Sensitive Drugs	Hydrogels protect from the UV light, enzymes, Ph changes make them perfect.
Biocompatibility and Mimicking of Biological Tissues	Biocompatible hydrogels copy natural tissue to encourage no immune response and increased biomedical acceptance by the body.
Minimal Invasiveness	Administration of hydrogels is non-invasive, Help for patient comfort and increased compliance with the treatment of eye, wound
Customisation and Flexibility	Formulating a complex polymer needs balancing stability, biocompatibility, and release characteristics.

**Table 1: Benefits<sup>(11)</sup>**

• **Limitations**

Key Features	Detailed Description & Therapeutics Impact
Complexity in Design:	Some cross-linkers, degradation products, and other stimuli-responsive elements may cause skin rash or allergic reaction.
Potential for Irritation and Allergic Reactions:	Some components of stimuli-responsive hydrogels' capacity can cause a rash on the skin or cause an allergy.
Limited Drug Loading Capacity:	Hydrogels' high moisture level prevents them from having a high dose of therapeutic agents available for drug delivery, and because networks within the hydrogel may be incompatible, additional therapeutic agents will not function effectively.
Challenges in Scaling Up Production:	Challenges related to making smart hydrogels at a large scale include developing appropriate methods

Difficulty in Regulatory Approval:	for quality control and consistency. Regulatory approval can be difficult for stimuli-responsive systems because they need real safety and efficacy testing, especially for novel polymer compositions and degradation products.
------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

**Table 2: Limitations<sup>(1)</sup>**

### Future Perspectives

Future hybrid hydrogels designed with AI-guided, 4D-bioprinted technologies through the use of organ-on-chip platforms will allow researchers to create personalised medicine. The capacity to take the place of the variable matrices characteristic of naturally derived materials (e.g., Matrigel) with standardised, programmable synthetic hydrogels can enable researchers to fine-tune their physical and biochemical signals, which supports the continued evolution of organoid technology by optimising the microenvironment for organoid cells.

### II. Conclusion

AI-guided engineering will help programmable hydrogels overcome the limitations of current manufacturing and drug-loading processes. The valuable synthetic microenvironments created by replacing variable matrices offer a way of connecting multifunctional materials with patient-specific clinical therapies.

### References

- [1]. Adina-Elena Segneanu, Ludovic Everard Bejenaru, Cornelia Bejenaru, Antonia Blendea, George Dan Mogosanu, Andrei Bita and Eugen Radu Boia. Advancements in Hydrogels: A Comprehensive Review of Nature and Synthetic Innovations for Biomedical Applications. *Polymers* 2025;17(15):2026.
- [2]. L. Gao, M. Ye, C. Wu. Cancer classification based on support vector machine optimized by particle swarm optimization and artificial bee colony. *Molecules*, 2017;22(12):2086
- [3]. B. Zhou, X. Li, Y. Pan, B. He, B. Gao. Artificial intelligence-assisted next-generation biomaterials: from design and preparation to medical applications. *Colloids Surf B Biointerfaces*, 2025: 255
- [4]. O. Ejeromedoghene, M. Kumi, E. Akor, Z. Zhang. The application of machine learning in 3D/4D printed stimuli-responsive hydrogels. *Adv Colloid Interface Sci*, 2025:336,
- [5]. W. Sheng, S. Cao, H. Qin, P. Liu, D. Wang, J. Weng, et al. Hydrogel in musculoskeletal disease: unravelling trends, research foci, and future trajectories via bibliometric insights (from 2000 to 2025). *Regen Ther*, 2025;30:544-557
- [6]. Guangli Xiang, Bohan Yin, Behzad Shiroud Heidari, George Youssef, Monika Gosecka, Mateusz Gosecki, Fernando G. Torres, Siu Hong Dexter Wong, Jagan Mohan Dodda. Programmable Hydrogels: Frontiers in Dynamic Closed-Loop Systems, Biomimetic Synergy, and Clinical Translation. *Advanced Science*. 2025;12(47)
- [7]. Qingquan Wang, Jiaqi Li, Nianping Feng, Yongtai Zhang. Smart responsive hydrogels combined with AI design for tumor therapy. *Acta Pharmaceutica Sinica B*. 2026
- [8]. S. Bhattacharya, K. Calar, P. de la Puente. Mimicking tumor hypoxia and tumor-immune interactions employing three-dimensional in vitro models. *J Exp Clin Cancer Res*. 2020;39(1):75
- [9]. L. Jiang, Y. Wufang, L. Xie, Y. Liu, X. Wu, et al. Experimental study on mechanism of stable drag reduction with hydrogel interface. *Tribol Int*. 2022;178
- [10]. Y. Liu, Y. Li, H. Shang, W. Zhong, Q. Wang, K. Mequanint, et al. Underwater instant adhesion mechanism of self-assembled amphiphilic hemostatic granular hydrogel from andrias davidianus skin secretion. *iScience*. 2022;25(10)
- [11]. Vibha Mishra, Dr. Abadhesh Kumar Niranjan, Vaishnavi Sahu. Smart Hydrogels: Mechanism, Recent Developments, and Potential in Drug Delivery System. *IJARESM*. 2024;12(80)